Doc no: CIBA/RA/Exp/004	Experiment				
Name of Department	Physics	_Location of Lab	S7-01-08		
Name of Laboratory	CIBA/ESP Nano Fabrication lab	_Name of PI	Asst Prof Jeroen Van Kan		
Name of Researcher/LO	Liu Fan & HueiMing	Name of Activity/Experiment	Poly-dimethylsiloxane (PDMS) casting/Bonding		

	1. Hazard Identification			2. Access the Risk			3. Risk Control			
No	Description/Details of Steps in Activity	Hazards	Possible Accident / III Health & Persons-at-Risk	Existing Risk Control (Mitigation)	Severity	Likelihood (Probability)		Additional Risk Control	Person Responsible	By (Date)
1	Stamp cleaning	IPA or acetone	IPA/acetone is irritating to respiratory system.	Wear mask and operating in fumehood	1	2	2			
2	Degassing PDMS prepolymer in desiccator	vacuum glass chamber	Heavy glass desiccators	Fixed glass desiccators on a secure surface	2	1	2			
3	Curing PDMS in oven	thermal oven	Thermal effects can burn skin	Setting thermometer to monitor the temperature	1	1	1			
4	Cleaning cured PDMS	IPA or ethanol	IPA is irritating to respiratory system	Wear mask and operate in fume hood	1	2	2			
5	Oxygen Plasma etching	UV light	Naked eye exposure	wear goggle when looking at the side chamber	1	1	1			

Conducted By HueiMing, Liu Fan, Ying Hui, John	Approved By	
	Name Asst Prof Jeroen Van Kan	
	Signature	
	First Revision Date: 31/10/2011	
	Approval date <u>31/10/2011</u>	Next Revision date 30/10/2014 (Maximum 3 years)